

Adiponectin with His Tag: Mouse, Recombinant Protein

Catalog Number:	1011-E4-ADPN
Product Specification:	Adiponectin with His Tag: Mouse Recombinant Protein
Species:	Mouse
Expression System	E. coli
Purity (by SDS-PAGE):	98%
Molecular Weight:	Da (calculated)
Size:	100 µg
Applications:	WB, ELISA, Cell culture
Formulation:	Prepared in PBS (1 mg/mL), Filtered (0.22 µm) & lyophilized.
Reconstitution:	Reconstitute at 1mg/mL in sterile water containing at least 0.1% human or bovine serum albumin.
Storage:	Store lyophilized and reconstituted proteins at -20°C for Long Term and at 4°C for < 2weeks. Avoid repeated freezing/thawing cycles.

Related Products:

1011-E1-ADPN	Adiponectin with <i>Flag</i> Tag: Human Recomb, in HEK293
1011-E2-ADPN	Adiponectin with <i>His</i> Tag: Human Recomb, in E.coli
1011-E3-ADPN	Adiponectin with <i>Flag</i> Tag: Mouse Recomb, in E.coli

Disclaimer: This product is for R&D use only, not for drug, household, or other uses. Dilutions and other information provided are recommendations only; the purchaser must determine the suitability of the product(s) for their particular use. Alamo[®] Laboratories Inc (Alamo[®] Labs) warrants that its product(s) conform to the information contained in this and other Alamo[®] Laboratories Inc publications. If the "Terms and conditions of Sale" are acceptable to you, then you may open the container(s) containing the product and, through such act of opening a container(s), will have shown your acceptance to "Terms and conditions of Sale". If the terms and conditions are not acceptable to you, then please return the container(s) unopened to Alamo[®] Labs for a complete refund of your payment.