

Mouse Skeletal Muscle Extract

Catalog Number:	2002-TEM-Mm
Product Specification:	Skeletal Muscle Extract
Species:	<i>Mouse</i>
Extraction Buffer:	Propriety
Format:	Supplied in 1X-Laemmli SDS-PAGE Buffer
Preservatives:	None
Size:	100 ul
Concentration:	1 mg/ml
Applications:	WB and ELISA Control Sample (Use in other applications, not tested)
Recommendations:	ELISA: Dilute 1:4 in 50 mM Carbonate Buffer pH 9.5 and use 100 ul/well as antigen. Western blot: Heat the mixture at 70 oC for 10 minutes and load ~10 ul/well on mini-gels.
Storage:	Store at -20°C for Long Term and at 4°C for < 1 week. Avoid repeated freezing/thawing cycles.

Disclaimer: This product is for R&D use only, not for drug, household, or other uses. Dilutions and other information provided are recommendations only; the purchaser must determine the suitability of the product(s) for their particular use. Alamo[®] Laboratories Inc (Alamo[®] Labs) warrants that its product(s) conform to the information contained in this and other Alamo[®] Laboratories Inc publications. If the “Terms and conditions of Sale” are acceptable to you, then you may open the container(s) containing the product and, through such act of opening a container(s), will have shown your acceptance to “Terms and conditions of Sale”. If the terms and conditions are not acceptable to you, then please return the container(s) unopened to Alamo[®] Labs for a complete refund of your payment.